



# FACE SHIELD

## Reusable facial protection for essential workers.



*Actual Face Shield may appear different than pictured.*

In March 2020, Actuated Medical, Inc. pivoted to support clinicians' needs for additional personal protective equipment (PPE) amidst a healthcare crisis. Our Face Shields are designed for multiple uses, allowing essential workers to disinfect and reuse.



Both Face Shield models use the same full coverage clear polycarbonate. Additionally, both allow room for glasses or goggles as well as additional PPE.

**Reinforced Headband**

**Simple Headband**

### Key Features

- + ANSI Z87.1 compliant
- + Clear polycarbonate for visibility
- + 300% thicker than typical face shields for reuse and durability
- + Open design reduces fogging
- + Chemical resistance to 1:8 bleach/water solution for over an hour
- + Disinfect according to your institution's protocols
- + Full-face protection
- + Multiuse for clinicians, police officers, security guards, cashiers, cooks, and more
- + One size fits most
- + Headbands may include Acrylonitrile Butadiene Styrene (ABS) and elastic

### Ordering Information

Minimum order: 20 Face Shields.

ITEM #	DESCRIPTION	QTY	PRICE / Shield
SFS-20	Reusable Face Shield with reinforced headband	20+	\$12 / ea. + s & h
EFS-20	Reusable Face Shield with simple headband	20+	\$6 / ea. + s & h

Net: 0.

**Call +1 (814) 355-0003 ext. 117**

**Email [FaceShield@ActuatedMedical.com](mailto:FaceShield@ActuatedMedical.com)**

These faceshields have not been FDA cleared or approved. The faceshields have been authorized by FDA under an EUA for use by healthcare providers as personal protective equipment, and are only authorized for the duration of the declaration that circumstances justifying the authorization of emergency use under Section 564(b)(1) of the Act, 21 USC 360bbb-3(b)(1) unless the authorization is terminated or revoked sooner. Actuated Medical, Inc. does not make any claims of safety or effectiveness of this product for any intended medical use.



Innovative motion + Positive outcomes

#### CONTACT SALES

Actuated Medical, Inc.  
+1 (814) 355-0003 x117  
[FaceShield@actuatedmedical.com](mailto:FaceShield@actuatedmedical.com)  
[www.ActuatedMedical.com](http://www.ActuatedMedical.com)



1100793368-003  
Pat. ActuatedMedical.com/ip