

# ACTUATEDMEDICAL® FACE SHIELD

Protecting all + Guarding futures

## Sustainably produced reusable face shields for additional coverage



Face shields help minimize respiratory droplet contact between people. They provide added eye protection while minimizing accidental face touching.

All Actuated Medical Face Shield models are manufactured using the same thick clear polycarbonate. The open top design reduces fogging but are available with an optional foam insert for over-the-top protection.

### Sustainable Protection

With face shields that have been demonstrated to last more than 5 months with daily use, **that's fewer face shields in the landfill**. Actuated Medical face shields are:

*“Here at Nestle we have over 1300 employees, some with bump caps or hard hats, and these face shields from Actuated Medical are the best option we’ve tried. They last through repeated uses and have better visibility than other brands.”*

+ Dennis Draleau, Occupational Health Nurse at Nestle Prepared Foods

#### Reusable



- ✓ Durable, Long-Lasting
- ✓ Easy to Clean

#### USA Strong



- ✓ Made in Pennsylvania
- ✓ USA-Sourced Materials

#### People Focused



- ✓ Full-Face Coverage
- ✓ Created New Jobs

**Longer lasting face shields = less plastic in landfills**

## Multiple Models, Same Quality Construction



Standard



Small



Hard Hat

*“These face shields are some of the best that I have ever used! They have very clear optics and do not fog up. I can even fit my dental loupes under them, which is critical for my work. They are easy to disinfect, which is particularly helpful right now, and are also longer and wider than normal shields, providing greater coverage. I got two to try them out initially but being so satisfied, I decided to order an additional ten shields.”*

+ Dr. Raya Ben Dor, Dentist

Actuated Medical Face Shields are made for adults and children alike. Specialty models for hard hats and select sports helmets are available. At Actuated Medical, we make shields that are:



### Reusable.

These face shields are constructed with 0.020” polycarbonate that has been proven durable for repeated use. Some customers have reported wearing the same face shield for more than 5 months at a time, reducing waste. They are easy to clean with simple dish soap and water or with a bleach/water solution.



### USA Strong.

Actuated Medical is a woman-owned small business located in beautiful Central Pennsylvania. Our plastic suppliers are also small businesses based in the United States, keeping manufacturing at home.



### People Focused.

Face shields provide splash protection for eyes and additional face coverage when paired with a cloth or medical face mask. Additionally, manufacturing face shields created additional jobs for Pennsylvania employees who were displaced by the pandemic.

## Key Features

- + Meets ANSI Z87 standard
- + Clear polycarbonate for visibility
- + 300% thicker than typical face shields for reuse and durability
- + Open design reduces fogging
- + Optional foam inserts close the open top for additional coverage
- + Chemical resistance to 1:8 bleach/water solution for over an hour
- + Easy to clean
- + Headbands may include Acrylonitrile Butadiene Styrene (ABS) and elastic

These face shields have not been FDA cleared or approved. The face shields have been authorized by FDA under an EUA for use by healthcare providers as personal protective equipment, and are only authorized for the duration of the declaration that circumstances justifying the authorization of emergency use under Section 564(b)(1) of the Act, 21 USC 360bbb-3(b)(1) unless the authorization is terminated or revoked sooner.

THE SHIELDS ARE PROVIDED “AS-IS,” AND ACTUATED MEDICAL DISCLAIMS ALL WARRANTIES, INCLUDING ANY (A) WARRANTY OF MERCHANTABILITY; OR (B) WARRANTY OF FITNESS FOR A PARTICULAR PURPOSE; WHETHER EXPRESS OR IMPLIED BY LAW, COURSE OF DEALING, COURSE OF PERFORMANCE, USAGE OF TRADE OR OTHERWISE. ACTUATED MEDICAL SPECIFICALLY DISCLAIMS THAT THE SHIELDS WILL PREVENT INFECTION OR THE TRANSMISSION OF VIRUSES OR DISEASES, INCLUDING COVID-19/CORONAVIRUS. ACTUATED MEDICAL FURTHER DISCLAIMS THAT THE SHIELDS WILL PREVENT PERSONAL INJURY OR DEATH. BUYER BEARS THE ENTIRE RISK AS TO THE PERFORMANCE, OPERATION, AND QUALITY OF THE SHIELDS.

Buyer acknowledges and understands that the shields have not been industry tested and may cause injury. Buyer agrees to indemnify, defend, and hold Actuated Medical harmless against any and all liability or loss arising in conjunction with or resulting from use of the shields (including any and all attorneys’ fees and expenses incurred by Actuated Medical).

**ORDER ONLINE:**

[www.FaceShield.ActuatedMedical.com](http://www.FaceShield.ActuatedMedical.com)

For larger quantities:

 **Call +1 (814) 355-0003 ext. 117**

 **Email [FaceShield@ActuatedMedical.com](mailto:FaceShield@ActuatedMedical.com)**



Innovative motion + Positive outcomes



1100793368-010  
Pat. ActuatedMedical.com/ip