


ACTUATED MEDICAL[®]
FACE SHIELD

Protecting all + Guarding futures


**NFHS compliant
 clear shield for most
 football helmets.**



Actual Face Shield may appear different than pictured.

Face coverings help minimize respiratory droplet contact between people. Compared to face masks, face shields offer an alternative covering for football players that provides clear visibility during practices and games.

Actuated Medical Football Helmet Face Shields feature the same 0.02 inch (0.5 mm) thick polycarbonate as its standard face shields for reusability and easy cleaning.



Top

Bottom

Attached to Helmet

Attached to Helmet

Key Features

- + NFHS approved
- + Full set contains 3 Top and 3 Bottom shields plus 30" Velcro[®]
- + Lowers-only set contains 3 Bottom shields plus 20" Velcro[®]
- + Clear polycarbonate for visibility
- + Universal design fits most types and sizes of football helmets
- + Chemical resistance to 1:8 bleach/water solution for over an hour
- + Replacement Velcro[®] available in multiple colors

Football helmet face shields include a top and bottom, with several attachment slots to fit most football helmet face masks. Velcro[®] must be cut to size (measurement included).

ORDER ONLINE:

www.FaceShield.ActuatedMedical.com

For larger quantities:

Call +1 (814) 355-0003 ext. 117

Email FaceShield@ActuatedMedical.com

These face shields have not been FDA cleared or approved. The face shields have been authorized by FDA under an EUA for use by healthcare providers as personal protective equipment, and are only authorized for the duration of the declaration that circumstances justifying the authorization of emergency use under Section 564(b)(1) of the Act, 21 USC 360bbb-3(b)(1) unless the authorization is terminated or revoked sooner.

THE SHIELDS ARE PROVIDED "AS-IS," AND ACTUATED MEDICAL DISCLAIMS ALL WARRANTIES, INCLUDING ANY (A) WARRANTY OF MERCHANTABILITY; OR (B) WARRANTY OF FITNESS FOR A PARTICULAR PURPOSE; WHETHER EXPRESS OR IMPLIED BY LAW, COURSE OF DEALING, COURSE OF PERFORMANCE, USAGE OF TRADE OR OTHERWISE. ACTUATED MEDICAL SPECIFICALLY DISCLAIMS THAT THE SHIELDS WILL PREVENT INFECTION OR THE TRANSMISSION OF VIRUSES OR DISEASES, INCLUDING COVID-19/CORONAVIRUS. ACTUATED MEDICAL FURTHER DISCLAIMS THAT THE SHIELDS WILL PREVENT PERSONAL INJURY OR DEATH. BUYER BEARS THE ENTIRE RISK AS TO THE PERFORMANCE, OPERATION, AND QUALITY OF THE SHIELDS.

Buyer acknowledges and understands that the shields have not been industry tested and may cause injury. Buyer agrees to indemnify, defend, and hold Actuated Medical harmless against any and all liability or loss arising in conjunction with or resulting from use of the shields (including any and all attorneys' fees and expenses incurred by Actuated Medical).



Innovative motion + Positive outcomes

